



ISOTHERMAL CONTAINERS

Code **BYF015 BYF016, BYF078**

Endo Plasticware isothermal containers, available in multiple capacities, have been developed to ensure the safe transport of temperature-sensitive substances while maintaining stable thermal conditions during use.

Features

- **Manufactured from high-quality polypropylene (PP)**, a material resistant to impacts, chemical agents, and temperature variations.
- **One-piece body with fully sealed bottom and walls**, ensuring a robust structure and excellent insulation capacity.
- **Bridge-style handles with safety locking system**, enabling secure transport and preventing accidental openings.
- **Exterior finished in blue with a white lid and interior**, allowing for easy identification and cleaning.
- **Operating temperature ranges from -40 °C to +80 °C**, suitable for demanding storage and transport conditions.
- **Not suitable for direct contact with food** – food products must be kept in their original packaging during storage and transport.



Technical Specifications

Code	Capacity (L)	Ext Dimensions (mm)	Weight (kg)
BYF015	25	420x250x350	2.0
BYF016	32	440x260x380	2.5
BYF078	45	520x310x400	4.2