

## Chemical waste disposal system, PE/PP

For collection of liquid chemicals in the laboratory. The inlet hopper made from PE-HD contains a self-closing float, overfill protection, and a splash guard. Additionally, a screw cap (GL 63) with sealing ring is included.

| Volume<br>l | Height<br>mm | Ø<br>mm | PU | Cat. No. |
|-------------|--------------|---------|----|----------|
| 10          | 560          | 222     | 1  | 151594   |



## Bowl, PP, with lid



White. Rectangular shape.

Broad, stable, easy to grip edge.

Especially easy to clean due to the rounded corners and edges and the smooth surfaces.

Suitable for contact with foodstuffs according to regulation (EU) No. 10/2011.

| Volume<br>l | L x W x H<br>mm | PU | Cat. No. |
|-------------|-----------------|----|----------|
| 17          | 430 x 331 x 195 | 1  | 43610    |



## Transport containers, PE-HD

Transparent.

Easy stackable.

With reinforcing ribs and integrated carrying handles.

Broad, stable edge.

| Volume<br>l | L x W x H<br>mm | PU | Cat. No. |
|-------------|-----------------|----|----------|
| 20          | 420 x 310 x 205 | 1  | 80602    |
| 46          | 600 x 365 x 260 | 1  | 80603    |
| 72          | 700 x 420 x 310 | 1  | 80604    |

