

## Adapter, Expansion

- Complies with DIN 12257 standard
- Large sized socket at the top and smaller sized cone at the base to connect dissimilar size ground glass joints
- Secure and airtight connection for distillation set ups

Item Code	Socket Size	Cone Size	Pack of
014.202.01	19/26	14/23	10
014.202.02*	24/29	14/23	10
014.202.03*	24/29	19/26	10
014.202.03A	29/32	14/23	10
014.202.04	29/32	19/26	10
014.202.05*	29/32	24/29	10
014.202.06*	34/35	19/26	10
014.202.07*	34/35	24/29	10
014.202.08*	34/35	29/32	10
014.202.09*	40/38	24/29	10
014.202.10*	45/40	29/32	10

## ASTM Sizes

Item Code	Socket Size	Cone Size	Pack of
014.470.01	19/22	14/20	1
014.470.02	24/40	14/20	1
014.470.03	24/40	19/22	1
014.470.05	29/42	24/40	1
014.470.09A	45/50	24/40	1

3.3  
BORO  
GLASS



1 x  
SOCKET

1 x  
CONE




ENQUIRE



ENQUIRE

\*Non DIN/ISO



## Adapter, Multiple joint with two parallel necks

- Complies with DIN 12594 standard
- Ground glass joints with a standard taper and smooth surface compatible with a variety of glassware to handle multiple liquids simultaneously

Item Code	Socket Size	Cone Size	Pack of
015.202.01	14/23	14/23	1
015.202.01A*	14/23	19/26	1
015.202.02	19/26	19/26	1
015.202.03*	19/26	24/29	1
015.202.04*	24/29	24/29	1
015.202.05*	19/26	34/35	1
015.202.06	29/32	29/32	1

## ASTM Sizes

Item Code	Socket Size	Cone Size	Pack of
015.470.01	14/20	14/20	1
015.470.03	19/22	19/22	1
015.470.04	24/40	24/40	1

3.3  
BORO  
GLASS



2 x  
SOCKET

1 x  
CONE



ENQUIRE



ENQUIRE

\*Non DIN/ISO

