

Volumetric Flasks, Amber, Class A

- Complies with DIN ISO 1042 Standards
- With Batch certificate
- With PE Stopper

Item Code	Capacity (ml)	N/S	Tolerance (±)	Pack of
130.202.01A	5	10/19	0.040	2
130.202.02AB	10	10/19	0.040	2
130.202.03A	20	10/19	0.040	2
130.202.04A	25	10/19	0.040	2
130.202.05A	50	12/21	0.060	2
130.202.05AB	50	14/23	0.100	2
130.202.06AB	100	12/21	0.100	2
130.202.06A	100	14/23	0.100	2
130.202.07A	200	14/23	0.150	2
130.202.08A	250	14/23	0.150	2
130.202.09A	500	19/26	0.250	2
130.202.10A	1000	24/29	0.400	2
130.202.11A	2000	29/32	0.600	1



ENQUIRE



3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

A
CLASS

PE



BATCH
CERTIFICATE

Trapezoidal Volumetric Flask, CLASS A, Amber Glass *New*

- Boro 3.3 glass
- With Batch certificate
- With PE Stopper

Item Code	Capacity (ml)	N/S	Pack of
135.202.00A	1	7/16	2
135.202.00BA	2	7/16	2
135.202.01A	5	10/19	2
135.202.02A	10	10/19	2
135.202.03A	20	10/19	2
135.202.04A	25	10/19	2
135.202.05A	50	12/21	2



ENQUIRE



3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

A
CLASS

PE



BATCH
CERTIFICATE