



## Pipeteurs

### Réf. - MBB002, MBB004, MBB007, MBG018

- The pipette controllers collection let you choose for the most effective solution in every application, helping you to get the most comfortable and fast manipulated of liquids through the use of graduated or volumetric pipettes.
- The pipette controllers as well as the pipette pumps, ergonomically designed, it have an internal mechanism that allows suck and dispense different volumes of liquids without change the manual position just with your thumb, avoiding the hand fatigue.
- With flexible connection and manufactured in rubber, they can be used with pipettes of different capacities.
- They are resistant to the action of acids and bases, and easily removable for cleaning.

## Manual controllers

### Réf. - MBB002, MBB004, MBB007

- Ergonomiques, légères et manieement facile avec une main
- Avec poire en gomme, manette pour aspiration et rejet du liquide et bouton pour rejet totale du liquide.
- Pour pipettes de 0.1 mL à 100 mL
- Avec filtre hydrophobe 3  $\mu\text{m}$  remplaçable pour protection du mécanisme interne



Référence	Couleur
MBB002	Bleu
MBB004	Vert
MBB007	Jaune

## Accessoires inclus

### Réf. - MBG018

Référence	Description
MBG018	Filtre pour pipeteur manuel DIGIPETTE