



**JIZERSKÁ PORCELÁNKA, s.r.o.**

Producer of laboratory and technical porcelains since 1947




## PRODUCT SPECIFICATIONS SHEET

### Plate dessicator porcelain

Jizerská porcelánka s.r.o. certifies that the specified products meet all specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

| Sign    | JIPO Code       | D (mm)  | Height (mm) |
|---------|-----------------|---------|-------------|
| 238/1   | 641 333 238 101 | 90 +3   | 5 +-1       |
| 238/3b  | 641 333 238 331 | 140 +-4 | 5 +-1       |
| 238/5a  | 641 333 238 521 | 189 +-4 | 5 +-1       |
| 238/8a  | 641 333 238 821 | 241 +-5 | 6 +-1       |
| 238/10a | 641 333 238 991 | 290 +-5 | 6 +-1       |



**Intended uses:** Laboratory applications

**Material:** Porcelain C110

**Shape:** Plate for dessicator

**Compliant to:** ISO material for the manufacture lies within the group C110 according to ČSN EN 60672 –3  
Product dimension tolerances are controlled by DIN 40680  
Acid resistance acc. to ISO 1775

**Properties:** Unglazed products can be used at temperatures of up to 1350°C / 2462°F and glazed ones at temperatures of up to 1000°C / 1,832°F  
Products have high resistance to acids and bases  
Density 2.4 g/cm<sup>3</sup>  
Dending strength of unglazed products (70MPa)  
Thermal conductivity (W.m-1.K-1) = 1.3



*JIZERSKÁ PORCELÁNKA, s.r.o.*

*Producer of laboratory and technical porcelainsince1947*



## PRODUCT SPECIFICATIONS SHEET

### Plate dessicator porcelain

---

The compliance of demands is proved by regular tests and followed noticed tests from an independent laboratory: Test according ČSN EN 60672-2 Determination of thermal shock resistance  $\Delta t$  230°C Test in liquid nitrogen Test resistance to acids and bases

**Disclaimer:** *Jizerská porcelánka s.r.o. wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use. If there's any doubt ask us.*