

## Short stem funnels

Code PNB001, PNB002, PNB003, PNB004, PNB005, PNB006, PNB007

Endo Glassware offers its range of smooth short stem funnels, designed for the safe transfer of liquids and filtration processes in laboratory environments. Manufactured from high-quality borosilicate glass, they provide excellent chemical and thermal resistance. Available in different sizes, they stand out for their durability and reliability in daily laboratory use.

### Characteristics

- Designed for the safe transfer of liquids and the filtration of solutions using filter paper
- Manufactured from borosilicate glass and in compliance with ISO 4798
- High chemical resistance to acids, bases, and most reagents commonly used in laboratories
- Excellent thermal shock resistance
- 60° cone angle, optimized to promote uniform flow and efficient filtration
- Smooth and homogeneous surface that facilitates cleaning, minimizes residue retention, and reduces the risk of cross-contamination



Code	Capacity (ml)	Funnel Ø (mm)	Stem Ø (mm)	Height (mm)	Recommended Filter Ø (mm)
PNB001	10	40	4	70	50-60
PNB002	30	50	7	95	50-60
PNB003	100	75	9	140	90-110
PNB004	150	90	10	170	110-125
PNB005	200	100	11	185	110-125
PNB006	400	120	12	220	125-150
PNB007	750	150	15	300	150-180