


Burette



Burette With PTFE Needle Valve Stopcock, Schellbach stripe, Class AS

- With Individual Certificate
- Complies with DIN/ISO 385 and USP standards
- Durable Blue Print
- Schellbach Stripe for Precise Liquid Dispensing and Accurate Meniscus Reading

Item Code	Capacity (ml)	Division	Tolerance ($\pm$ ml)	Pack of
112.223.01A	10	0.05	0.03	1
112.223.02A	25	0.10	0.05	1
112.223.03A	50	0.10	0.05	1
112.223.04A	100	0.20	0.10	1

3.3
BORO
GLASS

BLUE
PRINTING

Chemical
Resistant

NEEDLE

AS
CLASS

INDIVIDUAL
CERTIFICATE



ENQUIRE



Schellbach
Stripe



Burette with PTFE Needle Valve Stopcock, Schellbach stripe, Class AS

- With Batch Certificate
- Complies With DIN/ISO 385 and USP Standards
- Durable Blue Print
- Schellbach Stripe For Precise Liquid Dispensing and Accurate Meniscus Reading

Item Code	Capacity (ml)	Division	Tolerance ($\pm$ ml)	Pack of
112.202.01A	10	0.05	0.03	1
112.202.02A	25	0.10	0.05	1
112.202.03A	50	0.10	0.05	1
112.202.04A	100	0.20	0.10	1

3.3
BORO
GLASS

BLUE
PRINTING

Chemical
Resistant

NEEDLE

AS
CLASS

BATCH
CERTIFICATE



ENQUIRE