



# ENMIND EN-V5 Vet Infusion Pump

**Code ZMZ015**

AXAVET's infusion pump allows precise and safe administration of liquids and medicines. Quick access to the 20 most recent configurations.



- Small and lightweight
- Large touch screen for clear display of parameters including remaining infusion time
- Different settings available, such as: total perfusion volume, required flow rate, total perfusion time
- High sensitivity of settings (bubbles and occlusion pressure)
- High-performance battery (up to 10 hours of operation)
- Safety and operational alarms
- Self-diagnosis function
- Suitable for transfusions
- Compatible with standard droppers. Holder for serum carrier included

## Technical parameters

<b>Code</b>	ZMZ015
<b>Accumulated volume</b>	From 0,1 to 9999 mL
<b>Touch screen</b>	4,3 -inch color LCD
<b>Flow rate range</b>	0,1 – 2000 mL/h
<b>Modes:</b>	
Flow rate and volume	
Total time and volume	
Time and total flow rate	
Body weight	
Stored modes	
Micro mode	
<b>Accuracy</b>	±5%
<b>Air bubble detection</b>	20, 50, 100, 200, 400, 800 µl
<b>Occlusion settings</b>	5 levels (50, 150, 300, 600, 900 mm Hg)
<b>Data review</b>	5000 entries
<b>Protection class</b>	IP44
<b>Power supply</b>	100-240VAC, 50/60 Hz
<b>Dimensions (LxWxH)</b>	14x13x9 cm
<b>Weight</b>	1,45 kg