

Tubes

Digestion tube | Code PBV001, PBV002

- Designed for use in the determination of total nitrogen and phosphorus in automatic digesters.
- Made of borosilicate glass with thick, uniform walls, they are highly resistant to the action of strong acids and bases and to high temperatures.
- Round bottom



Code	Volume (mL)	Ø (mm)	Length (mm)	Constriction	For digester
PBV001	250	42	300	YES	Buchi
PBV002	250	42	300	NO	Tecator and others

Drying tube, ENDO GLASSWARE | Code PBV003



- Made of borosilicate glass

Code	Frosted
PBV003	14/23

Storage tube with hose connection, ENDO GLASSWARE | Code PBV004



- Made of borosilicate glass

Code	Capacity (mL)	Frosted	Ø hose connection (mm)	Total length (mm)
PBV004	20	14/23	10	127

Abductor tube, ENDO GLASSWARE

Code PBV005, PBV006, PBV007, PBV008, PBV009, PBV010, PBV011, PBV012, PBV023

- Made of borosilicate glass



Code	Ø (mm)	Angle (°)	Dimensions (mm)	Description
PBV005	8	90	100x100	-
PBV006	8	2x90	200x200x100	-
PBV007	7	90	100x100	-
PBV008	7	90	200x100	-
PBV009	7	2x90	200x100x100	-
PBV010	7	2x90	200x200x100	-
PBV011	7	2x90	100x210x200	With hook
PBV012	7	90,120	50x200x200	With hook
PBV023	7	2x90	200x210x100	With hook