



Plastic desiccators with button lid

Code **BDD020**, **BDD021**, **BDD022**

- Polypropylene (PP) base and disc and polycarbonate (PC) lid
- Lid with push button system facilitates hermetic seal
- Red coloured silicone gasket installed on the desiccator
- Transparent/translucent construction allows direct observation of the contents without the need to open the container



Technical specifications

Code	BDD020	BDD021	BDD022
Diameter	200 mm	250 mm	300 mm
Disc bores	3 mm	3 mm	15.5 mm
Dimensions	232x232x260 mm	276x276x310 mm	340x340x400 mm
Weight	0.85 kg	1.40 kg	2.40 kg