



Large capacity heating mantle with stirring

Code JGB009, JGB010

High-capacity equipment for heating liquids, mainly used in different industries, such as food, cosmetics, pharmaceuticals, petrochemicals, mining, etc.



Features

- Two models available: analogue (ref. JGB009) and digital (ref. JGB010)
- Agitation by Teflon-coated magnet
- Heating element: chrome-nickel wire
- Heating surface material: high-quality aluminium silicate fabric

Code	JGB009	JGB010
Capacity	20 L	
Maximum temperature	380°C	
Speed range	0-1600 rpm	
Heating power	2,5 kW	
Operating mode	Continuous	Continuous/Timed
Permissible ambient temperature	5-40°C	
Permissible ambient humidity	≤90%	
Power supply	200-240 V, 50/60 Hz	
Dimensions	Ø470x340 mm	600x457x360 mm
Approx. net weight	16 kg	