




Serological pipettes

Code **DGM020, DGM021, DGM022, DGM023, DGM024, DGM025**

- Made of medical grade polystyrene
- Gamma radiation sterilised
- Non-pyrogenic and DNase/RNase free
- Minimal adherence of the liquid to the pipette surface ensures precise delivery ($\pm 2\%$)
- With different coloured rings for easy identification
- Individually wrapped



Technical specifications

Reference	DGM020	DGM021	DGM022	DGM023	DGM024	DGM025
Volume (mL)	1	2	5	10	25	50
Length (mm)	270	278	340	340	311	354
External Ø (mm)	4,6	6	8,2	10	15	18
Graduation	1/100	1/100	0,1	0,1	0,2	1
Ring color	Yellow	Green	Blue	Orange	Purple	Black
Packing	100 pcs/bag	100 pcs/bag	50 pcs/bag	50 pcs/bag	25 pcs/bag	25 pcs/bag



Vol.
1 mL



Vol.
2 mL



Vol.
5 mL



Vol.
10 mL



Vol.
25 mL



Vol.
50 mL