



## BEAKER PP LOW FORM

Code **CBH001, CBH002, CBH003, CBH004, CBH005, CBH039, CBH040**

- Made of polypropylene.
- With spout to prevent spillage.
- With permanent embossed graduation.



CODE: CBH001



CODE: CBH003



CODE: CBH005

Code	Capacity	Graduation	Height	Diameter
CBH001	50 mL	10 mL	60 mm	42 mm
CBH002	100 mL	20 mL	70 mm	54 mm
CBH003	250 mL	50 mL	100 mm	71 mm
CBH004	500 mL	100 mL	120 mm	92 mm
CBH005	1000 mL	100 mL	145 mm	110 mm
CBH039	2000 mL	200 mL	185 mm	165 mm
CBH040	5000 mL	1000 mL	260 mm	215 mm