

Technical Specification

Reference: PCK401
Product: Pretreatment Pack

Description

PCK401 is a pretreatment pack including two cartridges: 20051 and 20331

Intended use:

- Feed water treatment and protection of reverse osmosis membrane
- Adsorption of impurities through activated carbon
- Removal of oxidation agents such as chlorine
- Particles retention
- Anti-scaling

Specifications



Ref	20051	20331
Use	Pretreatment filter	Pretreatment filter
Flow rate Max.	80 l/h	80 l/h
Application	Potable water	Potable water
Dimensions	Ø 74 mm – 340 mm	Ø 74 mm – 340 mm
Weight	290 g	600 g
Connections	¼" – ¼"	¼" – ¼"
Housing material	Polypropylene	Polypropylene
Contents	Sediment filter 5 µm	Granulated Active Carbon and anti scaling agent
Storage	Room temperature	Room temperature
Expiry date	Without expiry date	Without expiry date

This product is developed, produced and distributed according to management system that is certified for compliance with ISO9001.