



ISO BOTTLES ENDO GLASSWARE

CODE CFB001 to CFB004, CFB077 to CFB080

The ISO Endo Glassware bottles are manufactured in accordance with ISO standards and are suitable for the preparation, handling, and storage of laboratory solutions. Available in two types of glass, they provide options to meet different chemical and thermal resistance requirements.

General features:

- Blue cap and pouring ring made of PP, resistant to temperatures up to 140 °C
- Standard DIN GL45 thread
- White printed graduations and marking area
- Manufactured in clear glass according to ISO 4796



ISO ECO Bottles

- Manufactured from soda-lime glass

Code	Capacity (mL)	Base diameter (mm)	Height (mm)
CFB001	100	56	100
CFB002	250	70	143
CFB003	500	86	180
CFB004	1000	101	230



ISO Bottles

- Manufactured from borosilicate glass 3.3

Code	Capacity (mL)	Base diameter (mm)	Height (mm)
CFB077	100	56	100
CFB078	250	70	143
CFB079	500	86	180
CFB080	1000	101	230