



### ▶ Ultrasonic bath - Analogic

Ultrasonic cleaning baths allow all kind of instrument cleaning without detaching them. They are suitable for odontology and surgical material, jewellery, clockwork, optics and small tools, as well as for all kind of plastic, glass, porcelain, metal materials...

Cuenta con un cestillo de plástico donde alojar las piezas para su limpieza y tapa también en plástico, para evitar salpicaduras.



- 1| Plastic body and stainless steel tank.
- 2| Compact analogic ultrasonic water bath without heating for cleaning small items made of injected plastic with a stainless steel tank.
- 3| Provided with a plastic basket where placing the items to be cleaned and a plastic cover to avoid the splashing.
- 4| It includes an automatic disconnection device after 5 minutes of no use.
- 5| Especially intended for cleaning and disinfecting little instruments in dental offices, veterinary surgeries, hairdressers, etc.

Code	ZCC001
Mode	620/1
Capacity	0.6 L
Consumption	35 w
Useful dimensions	167x90x52 mm
External dimensions	223x114x123 mm
Weight	1 kg

### ▶ Cleaning solution



- 1| Easily mixed with water, it is formulated to remove general soils, fingerprints, stains and chemical residues from glassware and laboratory instruments.
- 2| Biodegradable, phosphate and caustic-free cleaner 20 mL per liter of water. Concentrated, divide it to 10-12%. Contains 1 L.

Code	ZCC004
Capacity	Contains 1 L