


Wire meshes / Wire triangles / Razors

WIRE MESH WITH FRAME

Stainless steel 1.4301, mesh: 0.6 x 2.0 mm,
fine mesh centre: heat resistant stainless steel


Edge length in mm	135	155	175	195
	80.247.935	80.247.955	80.247.975	80.247.995



WIRE MESH OUT OF STAINLESS STEEL WIRE

Stainless steel 1.4301, with folded up edge, mesh: 0.6 x 2.0 mm

Edge length in mm	115	135	155	175	195
	80.247.815	80.247.835	80.247.855	80.247.875	80.247.895



WIRE TRIANGLE

Stainless steel 1.4301, with clay pipes

Leg length in mm	40	50	60	70
	80.257.404	80.257.405	80.257.406	80.257.407



RAZOR

Long replaceable blade, incl. 10 spare blades

With plastic handle



80.066.071

RAZOR

Long replaceable blade, incl. 10 spare blades

With metal handle



80.066.072

SPARE BLADES

For 80.066.071 and 80.066.072

10 pieces

