

► Wide neck bottles PE

- 1 | Wide neck bottle, PE.
- 2 | Autoclavable up to 121°C.
- 3 | With cap.
- 4 | Natural colour.



CODE: CGD066



CODE: CGD068



CODE: CGD069

Referencia	Capacidad	Diámetro	Diámetro boca	Altura	GL
CGD066	50 ml	40 mm	24 mm	77 mm	32
CGD068	100 ml	48 mm	24 mm	92 mm	32
CGD067	250 ml	64 mm	31 mm	123 mm	40
CGD065	500 ml	75 mm	38.5 mm	153 mm	50
CGD069	1000 ml	94 mm	50 mm	202 mm	65
CGD070	2000 ml	120 mm	50 mm	243 mm	65