



WATER BATH

The water bath for education RS 3.5 has a cuve in polycarbonate with capacity of 3,5 litres, and with lacquered metal box.

- 1 | Heating by shielded heating element made of AISI 316 stainless steel
- 2 | Temperature regulation from ambient +5°C up to 90° C by bulb thermostat.
- 3 | Equipped with a safety thermostat that has a manual reset in accordance with DIN12877, class 2.
- 4 | Rack in stainless steel for 20 tubes Ø 16 mm.
- 5 | Illuminated ON switch.
- 6 | Indicator light showing heating element operation.



Code	JDB001
Model	RS 3.5
Volume	3.5 Litres
Temperature	RT +5°C > +90°C
Temperature regulation	Analogical
Work light	Yes
Cuve	Polycarbonate
Power	350 W