



Alcohol lamp

Code FFJ002, FFJ010, FFJ017

The Endo Glassware alcohol lamp consists of a rounded, flat-bottomed glass container with a screw-on lid and a cloth wick, which is in contact with the alcohol in the container. It is used to heat or sterilize samples or chemical reagents without requiring a high calorific value.

The thread has a safety system to prevent spillage in case of overturning.



Code	Volume (mL)	Total height (mm)	Ø wick (mm)	Wick length (mm)
FFJ002	100	90	4	120

Alcohol lamp wicks:



Code	Description
FFJ010	Wick metre Ø 4mm
FFJ017	Meter of cut wick Ø 4mm