



Specimen Container

Code ZMD010, ZMD011, ZMD012, ZMD013

- Made of high-quality polystyrene or polypropylene. Polypropylene containers can resist temperatures up to 121 °C and can be autoclaved
- Clear embossed graduation for easy reading
- Screw cap



Technical specifications

Reference	ZMD010	ZMD011	ZMD012	ZMD013
Material	PS	PP	PP	PP
Volume (mL)	30	60	90	120
Packaging	1000 pcs / box	500 pcs / box	400 pcs / box	250 pcs / box