

Extractors Soxhlet

Code PBC001, PBC002, PBC003

- The Soxhlet extractor is a glass apparatus used for the extraction of soluble components from solids using an organic solvent.
- Endo Glassware offers this series of Soxhlet extractors in various capacities and sizes, made of highly chemically and thermally resistant borosilicate glass, as part of their glassware equipment range.
- Each set includes a round-bottom flask, an extractor body, and a Dimroth condenser.



50 mL

Code	Description
PBC001	Extractor Soxhlet, 50 mL

Composed of:

Code	Description	Ground joint
PBN026	Dimroth condenser	29/32
PBR001	Extractor body	19/26 and 29/32
PBK003	Flat-bottom flask	19/26



125 mL

Code	Description
PBC002	Extractor Soxhlet, 125 mL

Composed of:

Code	Description	Ground joint
PBN027	Dimroth condenser	45/40
PBR002	Extractor body	29/32 and 45/40
PBK002	Flat-bottom flask, 250 mL	29/32



250 mL

Code	Description
PBC003	Extractor Soxhlet, 250 mL

Composed of:

Code	Description	Ground joint
PBN028	Dimroth condenser	55/44
PBR004	Extractor body	29/32 and 55/44
PBK010	Flat-bottom flask, 500 mL	29/32