




### SUB-SEMI SKIRTED, LOW PROFILE QPCR 96 WELL PLATE, 0,1 ML



- Ultra-thin walls to ensure maximum heat transmission
- Non-reactive inner surface and uniform wall thickness for consistency in error-free results
- Free of PCR inhibitors and endotoxins
- Certified DNAsa, RNAsa and human DNA free
- Designed in a special small size to avoid contamination between different investigations (10 plates per pack, x5 packs)

Code	Capacity	Wells	Steril	Pack
BGN028	0,1	96	No	10 uds x 5 packs

Dimensions of plates and wells. Microtube strip size.

