

## DATA SHEET

**Product:** Flat bottom tube PP, 50mL

**Reference:** BGD090



- Tube made of polypropylene (PP)
- Screw cap made of polyethylene (PE)
- With graduation (5mL, 10mL, 15mL, 20mL, 25mL, 30mL, 40mL, 50mL)
- External diameter: 30mm
- Total height (with cap): 108.5mm
- Temperature resistance: -80°C to 100°C
- Suitable for centrifugation. Max. speed: 9000 rpm    Max. RCF: 10000xg
- Presentation: Bag with 25 tubes