



▶ White rubber tubing

1 | Wide range of tubing for general uses in the laboratory; made of natural white rubber with a uniform color and total absence of bubbles or impurities.

2 | Slightly elastic and with a tension resistance up to 3 MPa, the white rubber presents a thermostability against temperature of boiling water that

avoid distortion or sticky texture after use.

3 | Among the main applications there are conduction of liquid or gas at low pressure and connection of Laboratory assemblies. They are supplied in rolls of 10 m.



Code	Inner Ø	Outer Ø	Thickness
FHN002	4 mm	6 mm	1 mm
FHN003	5 mm	7 mm	1 mm
FHN004	9 mm	12 mm	1.5 mm
FHN001	10 mm	14 mm	2 mm