



Chromatography Column With Sintered/Fritted Disc

- With PTFE needle Valve Stopcock
- Chemical Resistant With Porosity 0
- With socket for Connecting Reservoir Flasks
- Made from Boro 3.3 Glass

Item Code	Capacity (ml)	Body O.D. (mm)	Stem Length (mm)	Stem O.D. (mm)	Pack of
218.212.03	30	14	355	8	5

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant



ENQUIRE



Chromatography Column With Sintered/Fritted Disc

- With PTFE needle Valve Stopcock
- Chemical Resistant With Porosity 0
- With socket for Connecting Reservoir Flasks
- Made from Boro 3.3 Glass

Item Code	I.D. (mm)	Length (mm)	Socket Size	Pack of
220.202.01	10	100	14/23	1
220.202.02	10	200	14/23	1
220.202.03	10	300	14/23	1
220.202.04	10	400	14/23	1
220.202.05	20	200	19/26	1
220.202.06	20	300	19/26	1
220.202.07	20	400	19/26	1
220.202.08	20	500	19/26	1
220.202.09	26	400	24/29	1
220.202.10	26	500	24/29	1
220.202.11	30	400	29/32	1
220.202.12	30	500	29/32	1
220.202.13	30	400	34/35	1
220.202.14	30	500	34/35	1
220.202.15	40	400	40/38	1
220.202.16	40	500	40/38	1

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant

NEEDLE

1 x
SOCKET



ENQUIRE