

► Flacon PE à col large

- 1 | Col large, PE.
- 2 | Autoclavable à 121°C.
- 3 | Avec couvercle.
- 4 | Couleur naturelle.



RÉF: CGD066



RÉF: CGD068



RÉF: CGD069

Référence	Volume	Diamètre	diamètre col	Hauteur	GL
CGD066	50 ml	40 mm	24 mm	77 mm	32
CGD068	100 ml	48 mm	24 mm	92 mm	32
CGD067	250 ml	64 mm	31 mm	123 mm	40
CGD065	500 ml	75 mm	38,5 mm	153 mm	50
CGD069	1000 ml	94 mm	50 mm	202 mm	65
CGD070	2000 ml	120 mm	50 mm	243 mm	65