




► Evaporation dish with flat bottom, with spout



DIMENSIONS




- 1 | Made of borosilicate dish.
- 2 | They can be used in calcination processes, liquid evaporation or mixtures of acids and vapour.

Code	Capacity	d1.	d2.	h.
CNC001	15 mL	50 mm	25 mm	25 mm
CNC002	45 mL	60 mm	30 mm	30 mm
CNC003	90 mL	80 mm	35 mm	45 mm
CNC004	170 mL	95 mm	40 mm	55 mm
CNC005	320 mL	115 mm	50 mm	65 mm
CNC006	600 mL	140 mm	60 mm	80 mm

► Evaporation dish with round bottom, with spout



DIMENSIONS



- 1 | Made of borosilicate dish.
- 2 | They can be used in calcination processes, liquid evaporation or mixtures of acids and vapour.

Code	Capacity	d1.	d2.	h.
CNH001	35 mL	64 mm	25 mm	25 mm
CNC002	45 mL	60 mm	30 mm	30 mm
CNC003	90 mL	80 mm	35 mm	45 mm