



ENDO JARS

Code CCB001 to CCB010

- Made of polypropylene.
- With spout to avoid spillage and handle to ensure an easy grip.
- With permanent embossed graduation.
- Highly transparent, especially when containing liquid.
- Autoclavable

Graduated jug tall form, PP

Code CCB001 to CCB004



Code	Capacity mL	Graduation	d. mm	h. mm
CCB001	250	10/1	82	100
CCB002	500	10/1	95	125
CCB003	1000	100/1	110	180
CCB004	2000	20/1	155	220

Graduated jug short form, PP

Code CCB005 to CCB010



Code	Capacity mL	Graduation	d. mm	h. mm
CCB005	250	10/1	85	85
CCB006	500	10/1	90	118
CCB007	1000	50/1	119	130
CCB008	2000	200/1	143	167
CCB009	3000	100/1	167	180
CCB010	5000	250/1	190	270