

Tongs



CRUCIBLE TONGS

Stainless steel 1.4301, bent double, polished tips, electrolytically polished


Length in mm	200	220	250	300	400	500	600
	80.042.200	80.042.220	80.042.250	80.042.300	80.042.400	80.042.500	80.042.600



CRUCIBLE TONGS

Stainless steel 1.4301, bent double, tips unpolished

Length in mm	200	220	250	300
electrolytically polished	80.040.200	80.040.220	80.040.250	80.040.300
unpolished	80.840.200	80.840.220	80.840.250	80.840.300



CRUCIBLE TONGS

Stainless steel 1.4301, electrolytically polished, jointed with screw

Length in mm	200	220	250	300
bent double	80.081.200	80.081.220	80.081.250	80.081.300



BULB TONGS

Stainless steel 1.4301, electrolytically polished, cork inserts, for Ø 15 - 60 mm

Length in mm	230
	80.043.230



BEAKER TONGS

Stainless steel 1.4301, electrolytically polished, jaws wrapped with organic fleece, for Ø 60 - 150 mm,

Length in mm	300
	80.043.330