




Sterilized Bag Sealing Machine with roller holder

Code ZZC003

Equipment designed to provide a hermetic and secure seal on thermoplastic bags and materials. Its robust structure and advanced functionality ensure a perfect seal, protecting the contents from any leakage or contamination. Equipped with an integrated roll holder, it facilitates material handling, optimising the sealing process and guaranteeing consistent results with every use.

- Made of ABS combined with metal and with 4 rubber feet on the base
- Easy to operate, with pressure sealing and manual cut-off
- Start-up indicator light
- Preheat indicator light



* Roll not included

Code	ZZC003
Sealing width	15 mm
Max. sealing length	300 mm
Max. plastic thickness	0.4 mm
Dimensions	370x310x180 mm
Power	300 W
Weight	4.50 kg