

# Cylinder



## Mixing Cylinder with Hexagonal base, Class-A

- Complies with DIN EN ISO 4788 standards
- With conformity batch certificate
- With PE stopper

| Item Code    | Capacity (ml) | N/S   | Sub Division (ml) | Tolerance ( $\pm$ ml) | Pack of |
|--------------|---------------|-------|-------------------|-----------------------|---------|
| 142.202.00A* | 5             | 10/19 | 0.1               | 0.05                  | 2       |
| 142.202.01A  | 10            | 10/19 | 0.2               | 0.10                  | 2       |
| 142.202.02A  | 25            | 14/23 | 0.5               | 0.25                  | 2       |
| 142.202.03A  | 50            | 19/26 | 1.0               | 0.50                  | 2       |
| 142.202.04A  | 100           | 24/29 | 1.0               | 0.50                  | 2       |
| 142.202.05A  | 250           | 29/32 | 2.0               | 1.00                  | 2       |
| 142.202.06A  | 500           | 34/35 | 5.0               | 2.50                  | 2       |
| 142.202.07A  | 1000          | 45/40 | 10.0              | 5.00                  | 1       |
| 142.202.08A  | 2000          | 45/40 | 20.0              | 10.00                 | 1       |

3.3  
BORO  
GLASS



A  
CLASS



ENQUIRE

\* Non DIN



## Mixing Cylinder with Hexagonal base, Class-B

- Complies with DIN EN ISO 4788 standards
- With PE stopper

| Item Code    | Capacity (ml) | N/S   | Sub Division (ml) | Tolerance ( $\pm$ ml) | Pack of |
|--------------|---------------|-------|-------------------|-----------------------|---------|
| 142.222.00A* | 5             | 10/19 | 0.1               | 0.1                   | 2       |
| 142.222.01A  | 10            | 10/19 | 0.2               | 0.2                   | 2       |
| 142.222.02A  | 25            | 14/23 | 0.5               | 0.5                   | 2       |
| 142.222.03A  | 50            | 19/26 | 1.0               | 1.0                   | 2       |
| 142.222.04A  | 100           | 24/29 | 1.0               | 1.0                   | 2       |
| 142.222.05A  | 250           | 29/32 | 2.0               | 2.0                   | 2       |
| 142.222.06A  | 500           | 34/35 | 5.0               | 5.0                   | 2       |
| 142.222.07A  | 1000          | 45/40 | 10.0              | 10.0                  | 1       |
| 142.222.08A  | 2000          | 45/40 | 20.0              | 20.0                  | 1       |

3.3  
BORO  
GLASS



B  
CLASS



ENQUIRE

\* Non DIN