



Gloves

Code RBL002, RBL003

Silicone semi-glove

Code RBL002

- Material: Vapour phase silicone rubber
- Embossed inner surface
- Use: Handling of hot objects (glass, steel, porcelain, etc.)
- Maximum working temperature: 250°C (for prolonged exposure best below 200°C)
- Dimensions: 185x100 mm
- Specific weight: 1,04-1,05 g/ml
- Elongation at break: 300-400%
- Tensile strength: > 35.0 kg/cm²



Code	Size
RBL002	One size

Cotton glove

Code RBL003

- Material:
 - Backing: cotton (terry), cut/stitched.
 - Cuff: stretch fabric
- Cotton thickness: 850 g/m². Double knitted fabric.
- Use: Handling of hot objects of short duration.
- Maximum working temperature: 100°C



Code	Size
RBL003	One size