




Gauze iron wire

Code FFD001, FFD002, FFD003, FFD004, FFD005, FFD006, FFD007

Gauze iron wire with ceramic fibre core is used to optimise the heating process and ensure precise and safe results. They are manufactured by a combination of metallic materials and heat-resistant fibres, providing them with a number of advantages in terms of safety, stability and efficiency.

Technical characteristics:

- With ceramic fibre core.
- Asbestos-free.
- For heating substances with Bunsen burners, prevents breakage and direct contact with the flame.
- Even temperature distribution.



Code	Ø ceramic (mm)	Length (mm)	Width (mm)	Lot
FFD005	50	100	100	20
FFD001	80	125	125	20
FFD004	90	140	140	20
FFD002	100	150	150	20
FFD006	110	160	160	20
FFD007	130	180	180	20
FFD003	145	200	200	20