

## Adapter, Drying Tube

- Used to connect the drying/guard tube to a reaction vessel to keep the vessel free of moisture

| Item Code  | Cone Size | Pack of |
|------------|-----------|---------|
| 041.202.01 | 14/23     | 10      |
| 041.202.02 | 19/26     | 10      |
| 041.202.03 | 24/29     | 10      |
| 041.202.04 | 29/32     | 10      |

## ASTM Sizes

| Item Code  | Cone Size | Pack of |
|------------|-----------|---------|
| 041.471.01 | 14/20     | 1       |
| 041.471.03 | 24/40     | 1       |

3.3  
BORO  
GLASS



1 x  
CONE



ENQUIRE



ENQUIRE



## Adapter, Drying Tube Straight

- Used to connect the drying to a reaction vessel to keep the vessel free of moisture
- Standard joint for connecting the drying tube to a reaction flask

| Item Code  | Socket Size | Cone Size | Pack of |
|------------|-------------|-----------|---------|
| 041.470.02 | 19/22       | 19/22     | 1       |

3.3  
BORO  
GLASS



1 x  
CONE

1 x  
SOCKET



ENQUIRE

