

Adapter, Drying Tube

- Used to connect the drying/guard tube to a reaction vessel to keep the vessel free of moisture

Item Code	Cone Size	Pack of
041.202.01	14/23	10
041.202.02	19/26	10
041.202.03	24/29	10
041.202.04	29/32	10

ASTM Sizes

Item Code	Cone Size	Pack of
041.471.01	14/20	1
041.471.03	24/40	1

3.3
BORO
GLASS



1 x
CONE



ENQUIRE



ENQUIRE



Adapter, Drying Tube Straight

- Used to connect the drying to a reaction vessel to keep the vessel free of moisture
- Standard joint for connecting the drying tube to a reaction flask

Item Code	Socket Size	Cone Size	Pack of
041.470.02	19/22	19/22	1

3.3
BORO
GLASS



1 x
CONE

1 x
SOCKET



ENQUIRE

