

SPOON SPATULA WITH ROUND HANDLE

Stainless steel 1.4301



Length in mm	150	180	210	250
	80.210.150	80.210.180	80.210.210	80.210.250

MULTI-PURPOSE SPOON

Round handle Ø 5 mm,
spoon: approx. 15 x 35 mm, stainless steel 1.4301



Length in mm	150	180	210	250	300
for right hand	80.029.150	80.029.180	80.029.210	80.029.250	80.029.300
for left hand	80.329.150	80.329.180	80.329.210		

MULTI-PURPOSE SPOON

Round handle Ø 8 mm,
spoon: approx. 28 x 65 mm, stainless steel 1.4301



Length in mm	210	250	300	400	500
for right hand	80.029.218	80.029.258	80.029.308	80.029.408	80.029.508
for left hand	80.329.218		80.329.308		

SAMPLE SPOON

With filler pipe, stainless steel 1.4301



Length in mm	190
	80.029.190