




► Durham bell jar



- 1 | Our 1.5 ml Durham tube is made out of clear tubular glass of the 3rd hydrolytic class. This disposable glass tube has a round bottom.
- 2 | It's dimensions are $\varnothing 7.00 \times 50 \times 0.50$ mm.
- 3 | The Durham tubes are used in the microbiology to detect production of gas by microorganisms.
- 4 | These smaller test tubes are inserted upside down in another test tubes. At NAFVSM we offer the possibility to purchase the Durham tube already in another tube.

Code	BGG011
Type of glass	Sodalime glass
Type bottom	Redondo
Capacity	1.5 mL
Total length of product	50 mm
Wall tickness	0.50 mm
Diameter	7.00 mm