



ESPECTROSCOPIO DE MANO

HAND SPECTROSCOPE

REF : QLG008

INDEX OF LANGUAGE- ÍNDICE DE IDIOMAS

ESPAÑOL.....2
ENGLISH.....3

ESPECTROSCOPIO DE MANO

Sencillo instrumento para la examinación cualitativa de la comparación de espectros de absorción y emisión. Presenta los principales componentes como son colimador, prisma y telescopio montados dentro de un tubo recto, uno de cuyos extremos dispone de una rejilla ajustable para regular la entrada de luz y el otro extremo está provisto de un ocular. El extremo del ocular puede ser deslizado a lo largo de su longitud para enfocar el espectro. Suministrado completo en caja de madera.

HAND SPECTROSCOPE

A simple device for the rapid qualitative of the composition of emission and absorption spectra. It presents the main components as collimator, prism and telescope mounted inside a straight tube. One of the ends of the tube presents an adjustable slit to regulate the light entering inside and the other end is provided with an eyepiece. The drawtube of the eyepiece can be slide along its length for sharp focusing the spectrum.

Supplied complete in a wooden box.