

Noviglio, 03/07/2025

Item Reference:

Customer	Auxilab
Items codes	5571
Descriptions	URINE CONTAINER, STERILE SINGLE, 60 ml

Sterilization process:

Ray Emission	β -ray (Beta)
Minimum dose irradiated (kGy)	21,00.

Kartell S.P.A., Labware Division, declares under its own responsibility that the above mentioned products, to which this declaration refers, have been subjected to sterilization treatment. The process is validated UNI EN ISO 11137, UNI EN 11737, UNI EN 556, UNI EN 868, UNI EN ISO1167. It is guaranteed a SAL (Sterility Assurance Level) of 10^{-6}

Quality Assurance