




OLV-POC6

Portable oxygen concentrator



*Freedom to Breathe. Anytime. Anywhere.*

# Go Anywhere. Breathe Anytime.

Lightweight. Powerful. Reliable.

Designed for active lifestyles and real-life freedom.



## 4.4 lbs

**with Battery**

– Ultra-light for easy travel and daily use

## 6

**Pulse Levels**

– Equivalent to 6 L/min continuous flow

## 90%-96%

**Oxygen Purity**

– Stable, medical-grade oxygen output

## 9+

**Smart Alarms**

– Alerts for low oxygen, power, flow issues

## 1.5-5 hrs

**per Charge**

– Long-lasting battery for extended freedom

## Modular

**Design**

– Quick part swaps cut maintenance time & cost

# Smart, Portable, Ready for Real-Life Oxygen Needs

— Designed with Patients in Mind. Backed by Real-World Performance —



**Large LED display**  
— simple and clear operation



**Flowmeter Adjust**  
— Set oxygen level as needed



**On/Off Button**  
— one-touch operation



**Oxygen Outlet**  
— Secure and stable connection



**Logo**  
— Clean design,  
brand you can trust



**Inside Battery**  
— Backup power when  
you need it

Item	Details
Product Model	OLV-POC6
Flow Rate	6 pulse settings 1 fixed frequency setting
Oxygen Purity	90%-96% (1st-7th gears)
Dimensions	19.9×8.6×18.3cm (7.83×3.39×7.20 in)
Net Weight	1.98kg (4.37lbs) with battery
Operation Sound	≤60dB(A)
Warm-up Time	Approx. 2 minutes
Battery Operation Time	1.5-5H
Battery Charging Time	Within 2.5 Hours
Charging Voltage	AC100-240V50/60Hz DC19-21V
Warranty	One year

### Alarms & Safety Indicators

High Temperature Alarm	Low Pressure Alarm
Compressor Fault Alarm	No Breath Detected
Fan Malfunction Alarm	Warm-Up Reminder
Low Battery Warning	Power Failure Alarm
Low Oxygen Concentration Alarm	

### Standard Accessories

- 1 x Portable oxygen concentrator (main unit)
- 1 x AC power adapter
- 1 x Carrying bag
- 2 x Nasal oxygen cannulas
- 1 x HEPA filter





## **Zhengzhou Olive Electronic Technology Co.,Ltd.**

Add: 5th Floor, Block D., Building 18, Henan International University ,  
Science Park, (East District), 450001 Zhengzhou, Henan, China

Email: [Info@zzolive.com](mailto:Info@zzolive.com)

What's app: +86 18638288303

[www.zzolive.com](http://www.zzolive.com)

