


Issuer's name: **KAVALIERGLASS, a.s.**

Issuer's address: Křížová 1018/6, Prague 5
office: Sklářská 359, 285 06 Sázava, Czech Republic

Object of the declaration: **BOTTLE, SPECIFIC GRAVITY, GAY-LUSSAC**
Cat No. 1621/ 50 ml



Purpose of use: laboratory glassware

Material: **SIMAX®** Borosilicate glass

Eur. Ph. GLASS Type I

The object of the declaration described above is in conformity with the requirements of the following standards and regulations:

- ISO 3585 Borosilicate glass 3.3 – Properties.
- ISO 3507 Laboratory glassware -- Pyknometers
- Regulation EC No. 1935/2004 of 27 October 2004 on materials and articles intended to come into contact with food and repealing Directives 80/590/EEC and 89/109/EEC
- Regulation of Ministry of Health of Czech Republic No 38/2001, of the. Coll. of Laws, on hygienic demands for products intended to come into contact with food
- Regulation (EC) No 1907/2006 of the European Parliament and of the Council of 18 December 2006 - REACH - the object of the declaration does not contain any substance from the REACH - Candidate list SVHC
- Maximum permissible errors in capacity fulfill the values specified in Table 4/ISO 3507.

Table 4 - ISO 3507:

Nominal capacity, ml	1	2	5	10	25	50	100
Difference between actual capacity and nominal capacity, max., ml	±0,1	±0,3	±0,5	±1	±2	±3	±3

The producer confirms hereby that the characteristics, measures and accuracy of the products listed above are in full conformity with the provisions of the standard.

Additional information:

The producer declares that the products are safe when used in usual and proper way.

The producer has installed the Quality Assurance System according to ISO 9001 and thus guarantees that all products delivered to the market are in full conformity with the technical documentation and with all fundamental requirements to such products.

Certificate No. 04 100 940602 issued by TÜV CERT, Certification Body at TÜV NORD CERT GmbH.

Sázava, 14. 05. 2020

.....
Place and date of issue

Ing. Viera Košíková
.....
Quality manager