



# CERTIFICATE OF CONFORMITY

Issuer's Name / Producer : GLASSCO LABORATORY EQUIPMENTS PVT. LTD.

Issuer's Address / Producer : Manglai, P.O. Khudda Kalan, Ambala Cantt - 133004  
Haryana, India


Object of Declaration : Round Bottom Flasks

Material : Borosilicate Glass 3.3

Purpose of Use : Round-bottom flasks are used for uniform heating or boiling of fluid  
or chemical commonly used in laboratories.

## Specification of Round Bottom specified in Table-1, ISO 4797:2004

Catalogue No.	Product Description	Socket size	Capacity ml	Pack Qty.
057.202.04	Round bottom flask with joint	14/23	10	10
057.202.05	Round bottom flask with joint	14/23	25	10
057.202.08	Round bottom flask with joint	14/23	50	10
057.202.09	Round bottom flask with joint	19/26	50	10
057.202.10	Round bottom flask with joint	24/29	50	10
057.202.10A	Round bottom flask with joint	29/32	50	10
057.202.11	Round bottom flask with joint	14/23	100	10
057.202.12	Round bottom flask with joint	19/26	100	10
057.202.13	Round bottom flask with joint	24/29	100	10
057.202.14	Round bottom flask with joint	29/32	100	10
057.202.16*	Round bottom flask with joint	19/26	150	10
057.202.17*	Round bottom flask with joint	24/29	150	10
057.202.18*	Round bottom flask with joint	29/32	150	10
057.202.20*	Round bottom flask with joint	14/23	250	10
057.202.21	Round bottom flask with joint	19/26	250	10
057.202.22	Round bottom flask with joint	24/29	250	10
057.202.23	Round bottom flask with joint	29/32	250	10
057.202.26	Round bottom flask with joint	24/29	500	10
057.202.27	Round bottom flask with joint	29/32	500	10
057.202.28	Round bottom flask with joint	34/35	500	10
057.202.30	Round bottom flask with joint	24/29	1000	10
057.202.31	Round bottom flask with joint	29/32	1000	10
057.202.33	Round bottom flask with joint	24/29	2000	10
057.202.34	Round bottom flask with joint	29/32	2000	10



### \*Non ISO

The object of declaration described above is in conformity with the requirements of the following documents:

- DIN ISO 3585 Borosilicate Glass 3.3
- ISO 4797:2004 Laboratory Glassware - Round Bottom Flask

**Material:** The above are made from DIN ISO 3585:1998, heat resistant Boro 3.3 Glass that meet ASTM E-438, TYPE-1, CLASS- A Standard.

**Compliance:** These flask comply in every respect characteristics, measures and accuracy are in full Conformity With the provisions of the standard ISO 4797:2004 & IS 1381:2008.

### Additional Information:

Glassco has installed the quality control system as per ISO 9001 and products delivered to the market are in full conformity with the technical Specification as per standard and with all fundamental requirements to such products.

Certificate No. IPK 0039.06 issued by intertek certification limited, A UKAS Management System Accredited Body

This document has been produced electronically & is valid without signature

