



## AXL SEALERS

### Code ZZC004

- Digital model with display
- Made from a metal base with 4 rubber feet and a plastic top cover
- With a front-mounted handle for sealing and a side-mounted wheel for advanced the roll
- Integrated support tray



ZZC004

### Code ZZC005

- Made from a metal base with 4 rubber feet and a plastic top cover
- With front-mounted handle for sealing
- Indicator LEDs on the front panel
- Integrated support tray



ZZC005

### Code ZZC006

- Made entirely of metal with 4 rubber feet
- With a side-mounted handle for sealing
- Integrated support tray



ZZC006

### Code ZZC007

- Made from a metal base with 4 rubber feet and a plastic top cover
- With a side-mounted handle for sealing
- Without a support tray



ZZC007

Code	ZZC004	ZZC005	ZZC006	ZZC007
Model	SEAL-AXL 1	SEAL-AXL 2	Y400	S500
Max sealing length	300 mm	300 mm	250 mm	250 mm
Seal width	12 mm	12 mm	10 mm	10 mm
Display	Yes	No	No	No
Settings	Temperature: 80-240°C Sealing time: 2-8 s Sleep function: ON/OFF Counter reset function	Adjustable temperature thermostat	Adjustable temperature thermostat	Adjustable temperature thermostat
Power supply	220-240VAC 50/60 Hz	220-240VAC 50/60 Hz	220-240VAC 50/60 Hz	220-240VAC 50/60 Hz
Power	780 W	780 W	500 W	500 W
Dimensions	515x540x205 mm	470x540x205 mm	475x420x285 mm	420x310x260 mm
Weight	8.70 kg	7.60 kg	5.70 kg	4.75 kg