

DL-ZY04

Transfer electrophoresis tank



Technical Parameters

Effective transfer area: 110×90mm

Standard configuration

Name	Unit	Quantity
Lid with cables	Set	1
Bottom tank	pc	1
Transfer core	Set	1
Gel Holder Cassette	pc	4
Blue cooling unit	pc	1
Transfer Sponge Pad	pc	8

Features

- For fast, high-quality transfer of small protein condensation fronts
- The tank body is made of highly transparent polycarbonate material, which is one-time injection molding, no leakage and impact resistance
- Four 83*73mm coagulation fronts can be transferred at the same time
- The distance between positive and negative electrodes is 5.5cm, generating a strong electric field for efficient protein transfer
- 99.99% high-purity platinum electrode, better conductivity, corrosion resistance and aging
- Built-in blue cooling unit to quickly absorb the heat generated during electrophoresis
- Special open Gel Holder Cassette, easy to operate
- Equipped with a special bubble roller for easy removal of bubbles