



IRRIGATOR BAREL SHAPE WITH BENT STOPCOCK

Code **CLB002**, **CLB003** and **CLB004**

- Made of borosilicate glass
- Cover with knob and loose drop
- Tubulated near bottom



Code	CLB002	CLB003	CLB004
Ground	29/32	29/32	29/32
Volume	5 L	10 L	20 L
Ø lid (mm)	100 mm	127 mm	160 mm



LIDS FOR IRRIGATOR BARRELS

Code **CLB006** and **CLB017**

Code	Description
CLB006	Lid for irrigator barrel 5 L
CLB017	Lid for irrigator barrel 10 L