



Round basin - 4.5 L - white

Ref. : G606521



Plus point

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect.
Moulded rim makes it easy to grasp.

Product features

Customs code	39241000
EAN code	3100416065215
Net weight	0,36 kg
Height	95 mm
Diameter	320 mm
Capacity	4,5 L
Nesting rate	70,00 %
Product bottom internal height	90 mm
Product internal diameter	240 mm
Material	LDPE
Colour	White
French made	Yes
Without Bisphenol A	Yes

Logistic data

Number of pieces per pallet	450
Pallet weight	174 kg
Palette length	1200 mm
Pallet width	1000 mm
Pallet height	1630 mm
Number of units per layer	9
Number of layers per pallet	2
Number of pieces per sub-package	25
Weight in packaging	9 kg
Sub-package length	330 mm
Width under packaging	330 mm
Height under packaging	950 mm

All technical information on temperatures of use are given as an indication for the products concerned. The use of GILAC products remains the responsibility of the user.



Round basin - 7.5 L - white

Ref. : G606721



Plus point

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect.
Moulded rim makes it easy to grasp.

Product features

Customs code	39241000
EAN code	3100416067219
Net weight	0,46 kg
Height	125 mm
Diameter	350 mm
Capacity	7,5 L
Nesting rate	84,00 %
Product bottom internal height	120 mm
Product internal diameter	255 mm
Material	LDPE
Colour	White
French made	Yes
Without Bisphenol A	Yes

Logistic data

Number of pieces per pallet	540
Pallet weight	257,4 kg
Palette length	1200 mm
Pallet width	1000 mm
Pallet height	1340 mm
Number of units per layer	9
Number of layers per pallet	3
Number of pieces per sub-package	20
Weight in packaging	9,2 kg
Sub-package length	350 mm
Width under packaging	350 mm
Height under packaging	440 mm

All technical information on temperatures of use are given as an indication for the products concerned. The use of GILAC products remains the responsibility of the user.



Round basin - 9 L - white

Ref. : G606821



Plus point

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect.
Moulded rim makes it easy to grasp.

Product features

Customs code	39241000
EAN code	3100416068216
Net weight	0,65 kg
Height	125 mm
Diameter	390 mm
Capacity	9 L
Nesting rate	81,00 %
Product bottom internal height	120 mm
Product internal diameter	280 mm
Material	LDPE
Colour	White
French made	Yes
Without Bisphenol A	Yes

Logistic data

Number of pieces per pallet	240
Pallet weight	165 kg
Palette length	1200 mm
Pallet width	800 mm
Pallet height	1370 mm
Number of units per layer	6
Number of layers per pallet	4
Number of pieces per sub-package	10
Weight in packaging	6,5 kg
Sub-package length	390 mm
Width under packaging	392 mm
Height under packaging	330 mm

All technical information on temperatures of use are given as an indication for the products concerned. The use of GILAC products remains the responsibility of the user.



Round basin - 14 L - white

Ref. : G606921-1



Plus point

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect.
Moulded rim makes it easy to grasp.

Product features

Customs code	39241000
EAN code	3100416069213
Net weight	0,73 kg
Height	175 mm
Diameter	450 mm
Capacity	14 L
Nesting rate	89,00 %
Product bottom internal height	172 mm
Product internal diameter	270 mm
Material	LDPE
Colour	White
French made	Yes
Without Bisphenol A	Yes

Logistic data

Number of pieces per pallet	360
Pallet weight	274,8 kg
Palette length	1200 mm
Pallet width	1000 mm
Pallet height	1390 mm
Number of units per layer	6
Number of layers per pallet	3
Number of pieces per sub-package	20
Weight in packaging	14,6 kg
Sub-package length	445 mm
Width under packaging	415 mm
Height under packaging	400 mm

All technical information on temperatures of use are given as an indication for the products concerned. The use of GILAC products remains the responsibility of the user.