



Nitrile gloves

Code **RBL005, RBL006, RBL007, RBL008, RBL017**

- 100% unpowdered Nitrile glove. Highly elastic. Ambidextrous.
- For general use against minor hazards. Single use.
- Micro-textured outer surface on fingertips for better grip. Smooth inner surface to facilitate grip and reduce dermatitis.
- Reinforced cuff.
- Phthalate-free product

CE Certificate:

- MDR 2017/745 Class I. Related to Medical Devices
- EU 2016/425 Category III. Related to Personal Protective Equipment

Regulations regarding glove protection:

- EN ISO 21420:2020 General requirements for protective gloves
- EU 10/2021 Regulation on plastic materials and articles intended to meet foodstuffs, and its subsequent amendments
- EN 374 Protective gloves against chemicals and microorganisms (EN 374-1, EN 374-2, EN 374-4, EN 374-5 and EN 16523-1:2015)
- EN 455 Medical gloves for single use (EN 455-1, EN 455-2, EN 455-3, EN 455-4)
- ISO 13485:2016 International standard specific to medical devices
- ASTM D6978-05 Permeation test for chemotherapy drugs
- ASTM F1671 Resistance to viral penetration certification



Technical specifications

Code	Size	Weight	Colour
RBL005	S 6/7	5	Cobalt Blue
RBL006	M 7/8	5	Cobalt Blue
RBL007	L 8/9	5	Cobalt Blue
RBL008	XL 9/10	5	Cobalt Blue
RBL017	XXL 10/11	5	Cobalt Blue