

Medical / Face Shields / 701.00

701.00 - Clear Acetate

Code: MDU7001

Date: 15-12-25 1/2

701.00 face-shield's model provides protection against small flying particles characterized by an extremely low energy of impact while ensuring wearer's comfort. The thermoformed acetate visor will grant unobstructed vision (due to its anti-fog properties) even when user perspiration increases. The visor can be tilted up to an angle of 90° and can be replaced. This face-shield may be worn over ophthalmic spectacles and filtering half-masks depending on type and dimensions.



Frame mark (EU)

U EN166 S CE

Lens mark (EU)

U 1 S N CE

Frame technology

Overspec

Conditions of use

Standard brightness, Low brightness

Lens technology

Anti-Fog (N)

Lens color

Clear

Lens material

Acetate

Frame color

Grey/White

Standards

EN166

Mechanical resistance

Increased robustness (S)

Optional requirements

Resistance to fogging of oculars (N)

Optical class

Optical class (1)

Weight

68g

>> <https://www.univetsafety.com/en/products/medical/face-shields/701.00/701.00-clear-acetate/>



01.

01.

ANTI-FOG ACETATE VISOR



02.

02.

SOFT LATEX-FREE SPONGE ON THE
FOREHEAD

03.

ADJUSTABLE ELASTIC HEAD BAND

04.

VISOR INCLINABLE UP TO 90°

05.

INTERCHANGEABLE VISOR

06.

FACE SHIELD SUITABLE FOR USE IN
COMBINATION WITH PRESCRIPTION
GLASSES AND PROTECTIVE MASKS

Attention! Any images on this page illustrate the frame technologies of the product family in question and may represent frame and lens colours that do not correspond to the product code of the following sheet.

HEADQUARTERS

Univet S.r.l. Via Giovanni Prati 87 25086 Rezzato - Telephone 030 2499411 - Fax 030 2499430 - info@univetsafety.com
VAT. IT02506570163